

ASPIRATORE AD ANELLO UMIDO

PER ASPIRAZIONE CONTINUA SENZA NECESSITÀ DI ESSERE CONNESSI ALLA RETE IDRICA

WATER RING ASPIRATION

SEMI-DRY SUCTION UNIT: NO NEED TO BE CONNECTED TO THE WATER LINE

ASPIRATION À ANNEAU LIQUIDE

UNITE D'ASPIRATION AVEC SEPARATEUR AUTOMATIQUE: IL N'EST PAS NECESSAIRE D'ETRE CONNEXE A LA LIGNE D'EAU

BOMBA DE VACIO ANILLO DE AGUA

MÁQUINA DE ASPIRACIÓN Y UNIDAD DE SEPARACIÓN AUTOMÁTICA: NO HAY NECESIDAD DE ESTAR CONECTADO A LA RED DE AGUA

ASPIRATORI AD ANELLO LIQUIDO									Accessori				
	Kw	Volt	Hz	g/min	Vuoto: mbar	Portata lt/min	Dim		dB(A)	ITALIANO	ENGLISH	FRANCAIS	ESPAÑOL
RAIN 1	0,4	230	50	2800	150	1333	28 X 50 X 55cm 19 Kg	1	55	Centralina, antivibranti	control board, antivibrations	platine électronique, antivibrations	centralita, amortiguadores vibración
RAIN 1C	0,4	230	50	2800	150	1333	32 X 57 X 60cm 19 Kg	1	50	Centralina, antivibranti, carena insonorizzante	control board, antivibrations, metal cabinet	platine électronique, antivibration, capot metallique	centralita, amortiguadores vibración, carenado metálico
RAIN 2	0,94	230	50	2800	220	1500	34 X 43 X 54 cm 29 Kg	2	62	centralina, antivibranti	control board, antivibrations	platine électronique, antivibrations	centralita, amortiguadores vibración
RAIN 2+	0,6	230	50	2900 ~ 3720	200	1500	41 X 52 X 62 cm 31 Kg	1-2	65	Centralina, antivibranti, INVERTER	control board, antivibrations, INVERTER	platine électronique, antivibration, INVERTER	centralita, amortiguadores vibración, INVERTER
RAIN 3+	0,9	230	50	2900 ~ 3720	200	1700	41 X 52 X 62 cm 33 Kg	1-3	65	Centralina, antivibranti, INVERTER	control board, antivibrations, INVERTER	platine électronique, antivibration, INVERTER	centralita, amortiguadores vibración, INVERTER
RAIN 4+	1,1	230	50	2900 ~ 3720	200	1900	41 x 52 x 63 cm 37 Kg	1-4	68	Centralina, antivibranti, INVERTER	control board, antivibrations, INVERTER	platine électronique, antivibration, INVERTER	centralita, amortiguadores vibración, INVERTER



RAIN 1



RAIN 1C



RAIN 2



RAIN 2+ - RAIN3+ - RAIN 4+



4tek S.r.l
www.4-tek.it

V.le dell'Industria, 66
20037 Paderno Dugnano
(Milano) Italy
Ph: +39 02 91082211
Fax: +39 02 91084668

Email: info amministrative
riccardo@4-tek.it
info tecniche
stefano@4-tek.it
info commerciali
gabriella@4-tek.it